




STAR SCIENTIFIC GLASS CO.

(An ISO 9001-2000 Certified)

BEAKERS


S1000 BEAKERS

GRIFFIN, LOW FORM, WITH SPOUT

Capacity MI.	Approx OD mm.	Approx. Ht. mm	
5	20	25	
10	25	32	
25	34	45	
50	42	55	
100	50	73	
150	56	83	
250	68	95	
400	77	110	
500	83	115	
600	90	120	
1000	105	145	
2000	131	190	
3000	146	214	
5000	170	265	
10000	220	320	


S1040 BEAKERS

BERZELIUS, TALL FORM, WITHOUT SPOUT


Capacity MI.	Approx OD mm.	Approx. Ht. mm	
50	36	70	
100	50	80	
250	62	116	
500	75	137	
1000	90	190	

S1060 BEAKERS

BERZELIUS, TALL FORM, WITHOUT SPOUT

Capacity MI.	Approx OD mm.	Approx. Ht. mm	
100	50	80	
250	62	116	
500	75	137	
1000	90	190	

S1080 BEAKERS
PHILIPS, CONICAL, WITHOUT SPOUT

Capacity MI.	Approx OD mm.	Approx. Ht. mm	
250	75	123	
500	87	166	