




STAR SCIENTIFIC GLASS CO.


(An ISO 9001-2000 Certified)

BURETTES


S2118 BURETTES
FOR PINCHCOCK, WITH TIP CLASS B

Capacity MI.	Graduation Intervals ml	Tolerance \pm ml.	Height mm	
10	0.05	0.05	430	
25	0.1	0.1	470	
50	0.1	0.1	720	


S2121 BURETTES
ROTAFLO, FITTED WITH ROTAFLO GP (GENERAL PURPOSE) SCREW
THREAD
STOPCOCKS WITH PTFE KEYS, Class A

Capacity MI.	Graduation Intervals ml	Tolerance \pm ml	
10	0.05	0.02	
25	0.1	0.05	
50	0.1	0.05	
100	0.02	0.1	

S2122 BURETTES
 ROTAFLO, FITTED WITH ROTAFLO GP (GERERAL PURPOSE) SCREW
 THREAD
 STOPCOCKS WITH PTFE KEYS, Class B

Capacity MI.	Graduation Intervals ml	Tolerance \pm ml	
10	0.05	0.05	
25	0.1	0.1	
50	0.1	0.1	
100	0.2	0.2	

S2123 BURETTES
 STRAIGHT BORE STOPCOCK, CLASS A


Capacity MI.	Graduation Intervals ml	Tolerance \pm ml	Height mm	
1*	0.01	0.006	440	
+2*	0.02	0.01	440	
5*	0.02	0.01	590	
10	0.05	0.02	530	
25	0.1	0.05	530	
50	0.1	0.05	780	
100	0.2	0.1	830	

* with cup top

+ not covered in I.S.

S2124 BURETTES
STRAIGHT BORE STOPCOCK, CLASS B


Capacity MI.	Graduation Intervals ml	Tolerance ±ml	Height mm
1*	0.01	0.01	440
2*	0.02	0.02	440
5*	0.02	0.02	590
10	0.05	0.05	530
25	0.1	0.1	530
50	0.1	0.1	780
100	0.2	0.2	830



* with cup top

S212 BURETTES
DOUBLE OBLIQUE BORE STOPCOCK, 3 - WAY CLASS A

Capacity MI.	Graduation Intervals ml	Tolerance ±ml.	Height mm
10	0.05	0.02	530
25	0.1	0.05	530
50	0.1	0.05	780
100	0.2	0.1	830



S2126 BURETTES
DOUBLE OBLIQUE BORE STOPCOCK, 3 - WAY CLASS B

Capacity MI.	Graduation Intervals ml	Tolerance ±ml.	Height mm
10	0.05	0.05	530
25	0.1	0.1	530
50	0.1	0.1	780
100	0.2	0.2	830



S2127 BURETTES

DOUBLE OBLIQUE BORE ROTAFLO STOPCOCK 3 - WAY CLASS A

Capacity MI.	Graduation Intervals ml	Tolerance ±ml.	Height mm
10	0.05	0.02	530
25	0.1	0.05	530
50	0.1	0.05	780
100	0.2	0.1	830



S2128 BURETTES

DOUBLE OBLIQUE BORE ROTAFLO STOPCOCK 3 - WAY CLASS B


Capacity MI.	Graduation Intervals ml	Tolerance ±ml.	Height mm
10	0.05	0.05	530
25	0.1	0.1	530
50	0.1	0.1	780
100	0.2	0.2	830




S2129 BURETTES

STRAIGHT BORE PTFE KEY STOPCOCK CLASS A


Capacity MI.	Graduation Intervals ml	Tolerance ±ml.	Height mm
50	0.1	0.05	780




S2130 BURETTES
STRAIGHT BORE PTFE KEY STOPCOCK CLASS B

Capacity ML.	Graduation Intervals ml	Tolerance ±ml.	Height mm	
50	0.1	0.1	780	

S2145 BURETTES
DOUBLE OBLIQUE BORE ROTAFLO STOPCOCK 3 - WAY AUTOMATIC ZERO
CLASS A

Capacity ML.	Graduation Intervals ml	Tolerance ±ml.	Height mm	
10	0.05	0.02	530	
25	0.1	0.05	530	
50	0.1	0.05	780	


S2146 BURETTES
DOUBLE OBLIQUE BORE ROTAFLO STOPCOCK 3 - WAY AUTOMATIC ZERO
CLASS B

Capacity ML.	Graduation Intervals ml	Tolerance ±ml.	Height mm	
10	0.05	0.05	530	
25	0.1	0.1	530	
50	0.1	0.1	780	

S2147 BURETTES

DOUBLE OBLIQUE BORE STOPCOCK, 3 - WAY AUTOMATIC ZERO ACCURACY CLASS A.


Capacity MI.	Graduation Intervals ml	Tolerance \pm ml.	Height mm
10	0.05	0.02	530
25	0.1	0.05	530
50	0.1	0.05	780



S2148 BURETTES

DOUBLE OBLIQUE BORE STOPCOCK, 3 - WAY AUTOMATIC ZERO ACCURACY CLASS B.


Capacity MI.	Graduation Intervals ml	Tolerance \pm ml.	Height mm
10	0.05	0.05	530
25	0.1	0.1	530
50	0.1	0.1	780



S2149 BURETTES

AUTOMATIC ZERO, MOUNTED ON RESERVOIR CLASS A.


Capacity MI.	Graduation Intervals ml	Tolerance \pm ml.
10	0.05	0.02
25	0.1	0.05
50	0.1	0.05
100	0.2	0.1



Reservoir Capacity :

500 ml. for 10 ml. size and 2000 ml. for 25, 50 and 100 ml. sizes.

S2150 BURETTES
AUTOMATIC ZERO, MOUNTED ON RESERVOIR CLASS B.

Capacity MI.	Graduation Intervals ml	Tolerance \pmml.	
10	0.05	0.05	
25	0.1	0.1	
50	0.1	0.1	
100	0.2	0.2	

Reservoir Capacity :

500 ml. for 10 ml size and 200 ml for 25, 50 and 100 ml. sizes