



STAR SCIENTIFIC GLASS CO. BOTTLES


(An ISO 9001-2000 Certified)

S1220

BOTTLES

ASPIRATOR WITH OUTLET FOR TUBING

Capacity MI.	Approx OD mm.	Approx. Ht. mm	I.D. Neck mm
250	64	147	18
500	80	172	23
1000	104	200	28
2000	129	265	33




S1240

BOTTLES

ASPIRATION WITH OUTLET FOR STOPPER

Capacity MI.	Approx OD mm.	Approx. Ht. mm	I.D. Neck mm	I.D. Outlet mm
1000	104	200	28	23
2000	129	265	33	28
5000	190	320	44	28
10000	215	410	54	28




S1245

BOTTLES

ASPIRATOR WITH STOPPER AND STOPCOCK

Capacity MI.	Approx OD mm.	Approx. Ht. mm	I.D. Neck mm	I.D. Outlet mm
1000	104	200	B29	B24
2000	129	265	B34	B29
5000	190	320	B45	B29
10000	215	410	B55	B29
20000	300	485	B55	B29




S1250


BOTTLES

B. O. D. WITH STOPPER


Capacity MI.	Approx OD mm.	Approx. Ht. mm
125	55	120
300	70	150




S1260 BOTTLES
CENTRIFUGE, ORIGINAL FORM

Capacity MI.	Approx OD mm.	Approx. Ht. mm	I.D. Neck mm	
250	61	145	27	
500	73	182	27	


S1367 BOTTLES
MILK DILUTION, PLAIN, WITH SCREW CAP AND LINER

Capacity MI.	Approx OD mm.	Approx. Ht. mm	
160	150	19	

S136N BOTTLES
Mc CARTNEYFLAT WITH ALUMINUM SCREW CAP & RUBBER LINER, NEUTRAL

Capacity MI.	Approx OD mm.	Approx. Ht. mm	I.D. Neck mm	
210	64	40	150	


S1370 BOTTLES
MILK DILUTION, GRADUATED, SCREW CAP AND LINER

Capacity MI.	Approx OD mm.	Approx. Ht. mm	
160	150	19	

S1380
BAKELITE, WITH RUBBER LINER


SCREW

CAPS

Item Code	
S1380	


S1500
REAGENT, PLAIN NARROW MOUTH, WITH FLAT HEAD STOPPER

BOTTLES

Capacity MI.	Approx OD mm.	Approx. Ht. mm	Stopper	
60	45	93	B14	
125	55	120	B19	
250	64	147	B19	
500	80	172	B24	
1000	104	200	B29	
2000	129	265	B34	

S1501
REAGENT, WIDE MOUTH WITH SCREW CAP AND POURING RING

BOTTLES

Capacity MI.	Approx OD mm.	Approx. Ht. mm	I.D. Neck mm	
*30	30	75	30	
*60	45	90	30	
100	56	100	30	
250	70	138	30	
500	86	176	30	
1000	101	196	30	
2000	136	252	30	
5000	186	346	30	
10000	234	425	30	

* Without Pouring Rings


These Bottles are mechanically strong, chemically resistant and are provided with plastic pouring Ring for drip free operation. All bottles have only one size of screw Thread Caps and Pouring rings are Interchangeable and are made from P.P.

S1509

BOTTLES

REAGENT, AMBER, NARROW MOUTH WITH FLAT HEAD STOPPER

Capacity ML.	Approx OD mm.	Approx. Ht. mm	Size of stopper
60	45	93	B14
125	55	120	B19
250	64	147	B19
500	80	172	B24
1000	104	200	B29
2000	129	265	B34




S1585

BOTTLES

SOLUTION, PLAIN TOOLED NECK

Capacity ML.	Approx OD mm.	Approx. Ht. mm	I.D. Neck mm
250	64	147	18
500	80	172	23
1000	104	200	28
2000	129	265	33
3000	158	260	44
5000	190	320	44
9500	187	410	54
10000	215	476	54
13250	238	445	54
20000	300	485	54




S1589

BOTTLES

SOLUTION, AMBER, TOOLED NECK

Capacity ML.	Approx OD mm.	Approx. Ht. mm	I.D. Neck mm
250	64	147	18
500	80	172	23
1000	104	200	28
2000	129	265	33



S1625
RELATIVE DENSITY, WITH CAPILLARY BORE TEFLON STOPPER

BOTTLES


Capacity MI.	Approx OD mm.	Approx. Ht. mm	Stopper
10	24	46	B10
25	38	56	B10
50	49	66	B10



S1630
WEIGHING, WITH STOPPER

BOTTLES

Capacity MI.	Approx OD mm.	Approx. Ht. mm
5	20	40
15	25	60
*20	50	35
25	30	65
*40	60	40
60	40	90




* Squat Form

S1660
WASH, COMPLETE WITH STOPPER

BOTTLES


Capacity MI.	Stopper
250	B24
500	B29
1000	B29



S1661
ONLY FOR WASH BOTTLES OF 2771


STOPPERS

Capacity MI.	Stopper
250	B24
500	B24

1000	B29	
------	-----	---


S1662
ONLY FOR WASH BOTTLE OF 277

FLASKS

Capacity MI.	Approx. Ht. mm	Stopper	
250	135	B24	
500	175	B24	
1000	205	B29	


S1760
GAS WASHING, COMPLETE WITH STOPPER

BOTTLES

Capacity MI.	Stopper	
125	B29	
250	B29	
500	B29	

S1761
ONLY FOR GAS WASHING BOTTLES


STOPPERS

Capacity MI.	
125	
250	
500	

S1762
ONLY FOR GAS WASHING BOTTLES

BOTTLES

Capacity MI.	Approx OD mm.	Approx. Ht. mm	

125	50	155	
250	60	185	
500	70	245	