




STAR SCIENTIFIC GLASS CO.

(An ISO 9001-2000 Certified)

KETTLES

S6947 KETTLES RESIN REACTION


Capacity MI.	Approx OD mm.	Approx. Ht. mm
500	95	178
1000	108	178
2000	140	190
4000	140	343



The 500 ml and 1000 ml kettles have interchangeable covers with four openings to accommodate 24/29 interchangeable joints. Kettles of 2000 ml and 4000 ml capacity have interchangeable covers with openings to accommodate these 24/29 interchangeable joints and one 34/35 interchangeable joint.

S6948 BOTTOMS ONLY FOR 7058

Capacity MI.
500
1000
2000
4000



S6949 COVERS ONLY FOR 7058

Capacity MI.
500
100
2000
4000

